



B-200FLED Fluorescence Microscope



B-200FL/FL4 Fluorescence Microscope

- Wire winding device convenient for carry and storage (optional)



- Phase contrast unit, Independent phase contrast unit (apply to infinity optical system)



- Polarizer, analyzer for simple polarizing unit.



Analyzer installation

- Dry / Wet Dark Field Condenser



- Mirror



Polarizer installation

Simple polarizing unit
Polarizer, analyzer for simple polarizing unit

B-200F Fluorescence Microscope Outfits

Items	Specification	B-200F
Eyepiece	WF10× / 20mm	● ●
	WF10× / 20mm (Reticule 0.1mm,reticle)	○
	WF16×	○ ○
Infinity Plan Achromatic Objectives	PLAN 4× / 0.10	●
	PLAN 10× / 0.25	●
	PLAN 20× / 0.40(s)	○
	PLAN 40× / 0.66(s)	●
	PLAN 100× / 1.25(s, oil)	●
	PLAN 2.5× / 0.07	○
	PLAN 60× / 0.80(s)	○
Observing tube	Seidentopf Binocular Head	●
	Seidentopf Trinocular Head	○
Nosepiece	Inward Inclined positional quadruple nosepiece	●
Stage	Mechanical stage	●
Condenser	Abbe N.A.1.25	●
	Swing out achromatic condenser, N.A.0.9/0.13	○
Illumination device	6V/20W halogen (wide range input voltage: 100V~240V)	●
	6V/30W halogen (wide range input voltage: 100V~240V)	○
	3W-LED cold light illumination	○
Filter	blue (optional : yellow, green, gray)	●
C-mount	1×C-Mount (apply to digital camera)	○
	0.57×C-Mount (apply to camera system)	○
Adaptor for digital camera	Apply to CANON / NIKON / OLYMPUS etc.	○
Fluorescence Attachment	Epi fluorescence unit (B、G、Uv、V)	
	Epi fluorescence unit (B、G)	●
	100W mercury lamp house	●
	100Wmercury power source (optional: pointer or digital)	●
	Centering objective	●
	3W-LED Epi fluorescence unit (B、G)	○

Note: “●”In Table Is Standard Attachment. “○” Is Optional Accessories